

# Nephro Bag™ Nephro Cart™

NephroCan

SODIUM BICARBONATE CONTAINERS




# RETHINKING HEMODIALYSIS



2025

## CANADIAN LEADER IN HEMODIALYSIS



As a trusted Canadian provider, we are committed to delivering high-quality hemodialysis consumables, machinery, and equipment worldwide. Our goal is to ensure that our offerings are accessible and enhance patient care. We are dedicated to pioneering innovative therapies that revolutionize renal therapy, improving outcomes and enhancing the quality of life for patients globally.

The information and content provided in this brochure is the intellectual property of NephroCan and can only be used for purposes as intended by NephroCan. The content of this document shall not be altered, changed, or replicated. NephroCan claims possessory rights over the content of this document. The analysis and figures contained herein are based on NephroCan's interpretation and compilation of data and information of NephroCan's competitors in the hemodialysis market as published online and made available in the public domain.

# SODIUM BICARBONATE CARTRIDGES

Sodium Bicarbonate Cartridge Specifications	
Type	Sodium Bicarbonate Cartridge
Description	Sodium bicarbonate cartridges are designed ready to conduct on-line preparation of liquid bicarbonate solution for the safest and most effective hemodialysis treatment. External caps provide additional protection for hygienic disposal post-treatment.
Available Sizes	550g, 650g, 750g
Other	<ul style="list-style-type: none"> <li>» Sodium bicarbonate powder in a polypropylene cartridge for on-line production of liquid bicarbonate solution</li> <li>» Dust proof packaging</li> <li>» Open ports help eliminate perforation, reducing the risk of machine damage by unsolved powder</li> </ul>

## NEPHROCARD CHARACTERISTICS

Sodium bicarbonate cartridges, like the NephroCart™, play an essential role in hemodialysis by aiding in the correction of acid-base imbalances in patients with renal failure. These cartridges contain highly concentrated sodium bicarbonate in a dry powder form, encapsulated within dust-proof materials, and meet USP/BP and AAMI standards.

Available in three sizes, the NephroCart™ is universally compatible for on-line preparation of liquid bicarbonate, facilitating safe and efficient dialysis treatments. The cartridge is equipped with external caps to enhance protection and ensure hygienic disposal post-treatment, making the product a reliable choice for healthcare settings.



# SODIUM BICARBONATE BAGS

Sodium Bicarbonate Bag Specifications	
Type	Sodium Bicarbonate Bag
Description	Sodium bicarbonate bags are designed ready to conduct on-line preparation of liquid bicarbonate solution for one hemodialysis treatment, using significantly less storage space than sodium bicarbonate cartridges.
Available Sizes	550g, 650g, 750g
Other	<ul style="list-style-type: none"> <li>» Polyamide, polyethylene, hermetically sealed package ready for on-line preparation of liquid bicarbonate solution</li> <li>» Dust proof packaging</li> <li>» Angular design prevents loss of powder during use</li> </ul>

## NEPHROBAG CHARACTERISTICS

Sodium bicarbonate bags, such as the NephroBag™, are crucial in managing acid-base imbalances during hemodialysis for patients with renal failure. These bags feature highly concentrated sodium bicarbonate in a dry powder form, sealed within dust-proof materials to comply with USP/BP and AAMI standards.

Available in three sizes, the NephroBag™ is specifically designed for compatibility with a select range of hemodialysis machines, making the product suitable for all patients receiving treatment. The on-line preparation of liquid bicarbonate is seamless, ensuring safe and efficient dialysis treatments.



**NOTES:**

A series of horizontal dotted lines for taking notes, spanning the width of the page.

**NOTES:**

A series of horizontal dotted lines for taking notes, spanning the width of the page.



**NephroCan Inc.**

Suite 940-1040 West Georgia St. Vancouver, BC V6E 4H1  
Canada



**NephroCan S.r.l.**

Via Savoia 78 00198 Roma, Italy  
[www.NephroCan.com/EU-Rep](http://www.NephroCan.com/EU-Rep)



**Phone:**

+1 (604) 922-3110

**Fax:**

+1 (604) 922-3160

**Email:**

[info@NephroCan.com](mailto:info@NephroCan.com)  
[sales@NephroCan.com](mailto:sales@NephroCan.com)

**Website:**

[www.NephroCan.com](http://www.NephroCan.com)

